


AsepticPRO



No more operational failures during aseptic sampling:

With Norit Südmos fully automatic sampling system

- Sterile sampling with integrated CIP / SIP connection
- Single or multiple sampling possible
- Product viscosity up to 1,000 cP
- System pressure up to 10 barG
- Sterilisation temperature up to 150°C



System versions:

- Stand-a-lone Plug & Play with local control
- Integration of several systems with central control

Simple, secure and with an ergonomically safe operation!
Fit the sampling vessel - secure - and start sampling

Repeatable result...

...during continuous automatic process flow

Complete cleanability....

...through hygienic design and open flow cross section

Südmo Components GmbH
Industriestrasse 7 • D-73469 Riesbürg • Germany
T +49 9081 8030 • F +49 9081 803 158
E info@sudmo.de • I www.sudmo.de

Subject to technical revision 03/2008

E 1-2008

Norit Südmo Units

Norit
leading in purification

Südmo